


## Abstract

There is a large body of literature that documents the stresses that supraphysioxenic conditions place on human bone marrow mesenchymal stromal/stem cells (MSC), yet traditional room air culture practices with these clinically relevant cells still persist. We tested the hypothesis that MSC provided with physioxia (3 – 5% O<sub>2</sub>) just for incubation, and then handled in the supraphysioxia of room air (as in a biological safety cabinet), would have equivalent growth characteristics as MSC maintained in unbroken physioxia and CO<sub>2</sub>. Human bone marrow MSC cultured in triplicate T-75 flasks were housed for up to 8 passages in an external incubator fitted with an oxygen-controlled subchamber set to 5% O<sub>2</sub>/5% CO<sub>2</sub>. They were placed in HEPA-filtered room air conditions (Room-Air) for cell handling. The other set of cultures were housed within a closed processing chamber (Xvivo System) for full-time control of oxygen and CO<sub>2</sub>. Cell culture media were pre-equilibrated to the matching cell handling conditions before use. Cell growth was recorded using the CytoSMART cell imaging system. We found that the two sets of conditions did not produce equivalent cell growth. There were higher yields in each passage, and more total cell passages in cultures maintained full-time in physioxia conditions than in cultures handled under traditional supraphysioxenic conditions (two-tailed T test, unequal variances). We went on to look at how long it took pericellular and intracellular oxygen levels to recover after a quick room air medium change with medium that was not pre-equilibrated to physioxia levels. We measured pericellular oxygen levels with an oxygen probe and intracellular oxygen with an intracellular oxygen indicator dye. Pericellular oxygen levels in cultures handled in room air took over 80 minutes to equilibrate to the physioxia, far longer than it takes to modulate HIF-1 $\alpha$  levels in MSC. Intracellular oxygen took even longer to recover, over 150 minutes. Therefore, controlling oxygen levels around MSC during cell handling operations is critical to maximizing MSC yields both for research and clinical applications.


## Background

- Human mesenchymal stem/stromal cells (MSC) are of extremely high clinical value<sup>2</sup>.
- MSC originate from tissues with low O<sub>2</sub> levels and tend to use aerobic glycolysis<sup>1,5</sup>.
- HIF-1 $\alpha$  is affected within 5 minutes of O<sub>2</sub> change<sup>3</sup>
- MSC are negatively impacted by increased ROS in room air culture<sup>2</sup>.
- Supraphysiological and variable oxygen levels found in traditional cell culture and handling settings are highly stressful for MSC<sup>1</sup> and can affect proliferation<sup>2,4</sup>

## Methods




**Figure 1.** Lonza Poietics human bone marrow MSC cells were grown with Lonza hMSC Basal Medium plus MSCGM, hMSC Singlequot Kit, and Lonza Trypsin/EDTA for MSC, and Lonza DPBS w/o Ca/Mg (Lonza, Cologne, GmbH). Cells were thawed under physiologic conditions (37°C, 5%CO<sub>2</sub>, 5%O<sub>2</sub>) and immediately split into two conditions, A and B. One flask of each set was used for monitoring with the Lonza CytoSMART microscope and two flasks were used to pool with the first for cell counts at each passage. Cells were monitored for cell density and then subcultured when no more than 80% confluent. Medium was changed every 3-4 days even if cells did not need passaging. Cells were trypsinized, washed, and counted at each passage. Cell coverage data and time-lapse photomicrographs were provided by the CytoSMART. Cell counts were assessed for statistical overlap using Excel (Microsoft, USA).




**Figure 2.** Human bone marrow MSC (Lonza Gmb) were maintained in 3% O<sub>2</sub>, 5% CO<sub>2</sub> conditions for handling as well as incubation (A) in the Xvivo System. Incubators (black doors) open into a closed cell handling space with a heated floor and an atmosphere from tanked, filtered gases. To mimic the conditions in which cells are passaged in room air and returned to a tri-gas incubator, MSC were passaged at 20% O<sub>2</sub> plus 5%CO<sub>2</sub> and an immersion probe used to measure O<sub>2</sub> levels in the flask headspace and medium as the chamber was equilibrated to 3% O<sub>2</sub>. For intracellular O<sub>2</sub>, the MitoXpress Intra Intracellular Oxygen Assay (Agilent) was used. MSC grown at 3% O<sub>2</sub> were plated at three densities and loaded with the Intra reagent overnight. The next day, the cells were washed and the intracellular oxygen assays performed in the Xvivo System with an enclosed BMG plate reader.

## Results



**Figure 1.** Differential Effects of Full-Time Optimal and Part-Time Optimal Cell Culture Conditions on MSC. MSC grown under full-time optimal conditions for both incubation and handling consistently grew to higher yields than cells exposed to sub-optimal room air conditions for routine passage. This was consistent in early, mid-, and late-passage cultures, as assessed by CytoSMART cell coverage and by trypan blue cell counts. Cells in both conditions had good viability (inset), suggesting a cytostatic effect rather than a cytotoxic effect. Higher overall cumulative cell numbers were evident as well as more passages in unbroken optimal conditions. A 2-tailed paired T-test showed unlikely statistical overlap of yield at each passage (p=0.0079).



**Figure 2.** It Takes Almost Two Hours for a Vented T-75 Flask to Return to 3% O<sub>2</sub> After Passage in Uncontrolled Room Air Oxygen. Headspace and medium O<sub>2</sub> levels lagged far behind the chamber O<sub>2</sub> during equilibration in a T-75. A flask of cells returned to a tri-gas incubator does not equilibrate to chamber O<sub>2</sub> for hours. (Inset) Headspace and medium O<sub>2</sub> did not change when cells were handled under constant 3% O<sub>2</sub> and 5% CO<sub>2</sub>.

**Figure 3.** Cell Handling in Constant Conditions Eliminated Variability in Intracellular Oxygen. Intracellular O<sub>2</sub> levels were measured using Agilent's MitoXpress Intra Intracellular Oxygen Assay after a routine medium change at room air O<sub>2</sub> levels or at controlled 3% O<sub>2</sub>. Inside the cell, O<sub>2</sub> levels take hours to recover after an exposure to room air for routine medium changes. (Inset) Changing medium in a chamber with controlled O<sub>2</sub> and CO<sub>2</sub> conditions eliminated this variation.

## References

1. Boregowda, S. V., et al. (2012). Atmospheric oxygen inhibits growth and differentiation of marrow-derived mouse mesenchymal stem cells via a p53-dependent mechanism: Implications for long-term culture expansion. *Stem Cells*, 30(5), 975-987. doi:10.1002/stem.1069
2. Hoch, A.J., et al. Concise Review: Optimizing Expansion of Bone Marrow Mesenchymal Stem/Stromal Cells for Clinical Applications. *Stem Cells Translational Medicine*, 2(15), 4141. p. 412-412.
3. Jewell, U.R., et al. Induction of HIF-1 $\alpha$  in response to hypoxia is instantaneous. *FASEB J*, 2001, 15(7), p. 1312-4.
4. Paquet, J., et al. (2015). Oxygen Tension Regulates Human Mesenchymal Stem Cell Paracrine Functions. *Stem Cells Transl Med*, 4(7), 899-921. doi:10.5966/sctm.2014-0190
5. Phinney, D.G., et al. Mesenchymal stem cells use extracellular vesicles to outsource mitogenesis and shuttle microRNAs. *Nat Commun*, 2015, 6, p. 8472.

## Conclusions

- Cell handling in room air O<sub>2</sub> produces out-of-optimum pericellular and intracellular conditions for hours after cells are back in the incubator. Cell handling in full-time controlled conditions prevents these variations.
- Full-time control of oxygen for incubation and cell handling increases MSC yields over part-time optimal conditions.