

# FULL-TIME PHYSIOXIC CULTURE CONDITIONS PROMOTE MSC PROLIFERATION AT PHYSIOXIC CONDITIONS BETTER THAN HYPOXIC PRE-CONDITIONING

Alicia D. Henn, Yan He, Shannon Darou, Randy Yerden  
BioSpherix, Ltd.

## ABSTRACT

Cells that are cultured and handled in traditional hyperoxic room air experience a dramatic shift in  $O_2$  levels when injected *in vivo*. It has been reported that hypoxic pre-conditioning of mesenchymal stromal cells (MSC) in culture can increase cell survival *in vivo* after administration. Hypoxic pre-treatments subject the cells to lower  $O_2$  atmospheric conditions for lengths of time that range from 15 minutes to 36 hours, with another period of high oxygenation in room air conditions before administration. We previously showed that Human bone marrow MSC yields are improved by enclosing cell-handling processes as well as incubation for a full-time controlled physioxic atmosphere. Here, we used the Xvivo System (BioSpherix), with fully-enclosed and environmentally-controlled cell incubation and handling chambers, to compare different hypoxic pre-conditioning regimens to full-time physioxic conditions. Our null hypothesis was that the different culture conditions make no difference to MSC proliferation rates. MSC were incubated for 12-36 hours at 1%  $O_2$  before return to traditional supraphysioxenic room air incubator oxygen levels for 24 hours. The cells were then incubated at 5%  $O_2$  (venous levels) as if injected *in vivo*. These conditions were compared with unbroken 5%  $O_2$  for both cell handling and culture. Using an immersion  $O_2$  probe, we recorded vessel headspace and pericellular medium  $O_2$  levels during each regimen. We looked for changes in cell division rates. MSC pericellular  $O_2$  levels lagged far behind atmospheric  $O_2$  level changes, so MSC experienced low  $O_2$  conditions for far less time than the pre-conditioning period. Unbroken  $O_2$  conditions were more favorable for keeping MSC in the cell cycle than any pre-conditioning regimen, disproving the null hypothesis. We concluded that hypoxic pre-conditioning is not as favorable for MSC yields as full-time physioxic conditions for cell handling and incubation.

## BACKGROUND

- Exposing MSC cultures grown in room air to low  $O_2$  has been shown to benefit MSC function<sup>1</sup>
- Low  $O_2$  induces HIF-1 $\alpha$  stabilization and activation of downstream genes involved in cell proliferation<sup>2</sup>
- We previously showed that human bone marrow MSC proliferate for more passages at full-time physioxic conditions than in room air culture<sup>3</sup>


## OBJECTIVES

- Assay MSC proliferation rates after different times at 1%  $O_2$  for hypoxic preconditioning and compare with full-time physioxia (5%  $O_2$ )
- Assay final cell numbers after treatment
- Assess pericellular oxygen levels at each condition

## REFERENCES

1. Bader AM, Klose K, Bieback K, Korinth D, Schneider M, Seifert M, Choi Y-H, Kurtz A, Falk V, Stamm C: Hypoxic preconditioning increases survival and pro-angiogenic capacity of human cord blood mesenchymal stromal cells *in vitro*. *PLoS one* 2015, 10(9):e0138477.
2. Lambertini E, Penolazzi L, Angelozzi M, Bergamin LS, Manfredini C, Viecei Dalla Sega F, Paoletti F, Lisignoli G, Piva R: Hypoxia Preconditioning of Human MSCs: a Direct Evidence of HIF-1 $\alpha$  and Collagen Type XV Correlation. *Cellular Physiology and Biochemistry* 2018, 51(5):2237-2249.
3. Henn A, Farrell G, Darou S, Yerden R: Shock from Exposure to Room Air Conditions Alters Individual Mesenchymal Stem Cell Fate, Population Dynamics, and Batch Yields. *International Society of Cell Therapy North America* 2017.




## EXPERIMENTAL DESIGN



**Figure 1. Experimental Set-Up.** (A) Schematic, dynamic  $O_2$  control for cell culturing and handling. Human bone marrow MSC in T-25 flasks were cultured under 5% constant  $O_2$  (physioxia; black line), 18% constant  $O_2$  (blue line), or 12, 24 or 36 hr of 1%  $O_2$  (hypoxic preconditioning; orange solid line). The cells were harvested and plated into a 96-well plate and cultured at 5% constant  $O_2$  as if injected *in vivo*. (B) MSC were cultured, handled and observed in fully controlled conditions at all times using the Xvivo System. Each process chamber or incubator (black doors) had individually controlled temperature,  $O_2$ ,  $CO_2$  and relative humidity (RH).



## RESULTS



## CONCLUSIONS

**Pericellular Oxygen Conditions Lag Behind Incubator Atmospheric Oxygen Changes**  
**Hypoxic Preconditioning Increased Cell Proliferation as Compared with Room Air Conditions, but Not as Compared to Full-Time Physioxic Conditions**  
**Overall Cell Growth is Negatively Affected by Room Air Conditions, so Physioxic Conditions are Better**